

Use case scenario Process

Nr.	Wer	Was
1.0	Patient	Patient confirms pressing the Display
1.1	Display	Sends confirmation to the PIP
1.2	PIP	PIP sends command to quit alarm
1.3	Alarm	Alarm quits and sends confirmation back to PIP
1.4	PIP	Alarm quits and sends confirmation back to PIP
1.5	Display	Battery-Live will be calculated and information sent to the display
1.6	Display	Display shows remaining battery-Livetime
1.7	Patient	Confirms action
1.8	Display	PIP receives Patient's confirmation