Use case scenario

Name	Batterie-Alarm
Nummer	1
Kurzbeschreibung	Patient kann den Batterie-Alarm mittels Druck auf das Display D1 quittieren.
Akteure	Patient, Display, Alarm, PIP
Auslöser / Vorbedingung	Alarm ist Ausgelöst
Ergebnis / Nachbedingung	Alarm ist Ausgeschaltet


Use case scenario Process

Nr.	Wer	Was	
1.0	Patient	Patient confirms pressing the Display	
1.1	Display	Sends confirmation to the PIP	
1.2	PIP	PIP sends command to quit alarm	
1.3	Alarm	Alarm quits and sends confirmation back to PIP	
1.4	PIP	Alarm quits and sends confirmation back to PIP	
1.5	Display	Battery-Live will be calculated and information sent to the display	
1.6	Display	Display shows remaining battery-Livetime	
1.7	Patient	Confirms action	
1.8	Display	PIP receives Patient's confirmation	

Alternative Use case scenario

Nr.	Wer	Was
1.7.1	Patient	Doesn't accept action
1.7.2	Display	Display failure message
1.7.3	PIP	Send Reminder after α-sec
1.7.4	Display	Displays reminders
1.7.5	Patient	Patient confirms action else back to 1.7.1

Acitivity Diagram Elaboration



Acitivity Diagram Elaboration

